

PROMOTIONAL OFFERS*

Valid from 3th of February till 30th of June 2025

Only for participating dealers Benelux.
All prices are in €, exclusive VAT.



*Benelux area



CONTACTS

Area Sales Manager Belux Equipment & Pharma

Raphaël NEF

Tel.: +32 (0) 472 711 580
raphael.nef@acteongroup.com

Area Sales Manager Imaging The Netherlands

Andy LUBERTI

Tel.: +31 (0) 613 948 984
andy.luberti@acteongroup.com

Administration Commerciale

Sabine PARIS

Belgium / Nordics
Tel.: +33 (0) 556 349 332
Fax: +33 (0) 566 349 331
sabin.paris@acteongroup.com

Nirina INISAN

Luxembourg / The Netherlands
Tel.: +33 (0) 556 349 250
Fax: +33 (0) 566 349 331
nirina.inisan@acteongroup.com

ACTEON® CUSTOMER SERVICE

ACTEON's* customer service
can be contacted at:

Imaging

PSPIX® - SOPIX® - CAMERAS - X-MIND®

Tel.: +33 (0) 481680781
cs@acteongroup.com

X-MIND® trium - X-MIND® prime

Tel.: +39 (0) 331 37 67 62
imaging.italysupport@acteongroup.com

Equipment / Pharma

Tel.: +33 (0) 556 349 163
cs@acteongroup.com

ACTEON®

17 avenue Gustave Eiffel
33700 MERIGNAC - FRANCE
TEL: +33 (0) 556 340 607
FAX: +33 (0) 556 349 292
info@acteongroup.com
www.acteongroup.com

Non contractual document
Copyright© 2025 ACTEON®. All rights reserved. No information or
part of this document may be reproduced or disclosed in any format
without the prior authorisation of ACTEON®.
For professional dental use only.
Read the instructions in the notice carefully before use.
All information that is essential to ensure correct use of these
devices is included in the user manual or the summary of the
product characteristics available on www.acteongroup.com.
Creation date: January 2025

Our Values

A human focused, service-oriented company

The customer knows best: ACTEON®'s solutions are inspired and guided by the specific need of its diverse partners, be they healthcare professionals or distributors. ACTEON®'s long-standing presence in the market has enabled the company to develop granular understanding of each type of customer and to create best-fit solutions.

Innovator at heart

Creative thinking has always been a part of ACTEON®'s DNA and is still a core pillar of the company today. ACTEON® is building on this strength everyday to set new standards for the industry and improve people's lives.

Excellence and expertise

ACTEON® has always stood for high quality and high performance. From advanced technologies to the manufacturing of manual instrumentation, the company has made a name for itself and grown trust through its requirements and its extensive experience in the industry.

Entrepreneurial spirit


ACTEON® has created a work environment that enables collaborators to make the most of their capabilities, empowers them and encourages creativity and action.

Cooperation & knowledge mindset

Team spirit is playing a central role in connecting the different entities of the group. It also infuses ACTEON®'s approach to its customers and partners and is a guiding principle in its ambition to transform the industry and everyday practice.

Exclusive and patented Color Coding System™

Power selection is accurately determined with the Color Coding System™, which intuitively associates each tip color ring with one of the 4 power ranges



newtron p5XS


1
Ref. F00246



1S
Ref. F00245



H3
Ref. F00369



P5XS b-LED incl. with the following tips:
Ref. F62100

~~€2788,25~~
€2099




N°4
Ø22
Ref. MP3622PH

NOW INCLUDED

PURE REFLECT®
with a value
of €54,89

Ref. MP3622PH

Newtron®



Linear and regular vibrations

- Painless treatment for the patient



Pièce à main fine et légère

- Decreased fatigue for the clinician




Pas de chaleur ni de vibration de la pièce à main

- Better tactile sense



Perfectly controlled vibrations

- Precise & gentle treatments




Specific steel quality

- Anatomy & enamel preservation




Total irrigation control

- No heat is safer for tissues




Automatic and continuous frequency adjustments

- Maximal and constant efficiency whatever the load applied



Real-time power compensation

- Continuous performances
- No force needed



Powerful cavitation

- Deposits fragmentation
- Disinfecting effect

Newtron® handpiece



ref. F12281 - Autoclavable

~~€635,50~~

€339

Newtron® LED handpiece



ref. F12609

~~€879,63~~

€495

Newtron® Slim B. LED handpiece



- with blue LED ring: ref. F12900
- with white LED ring: ref. F12905

~~€879,63~~

€549

Hygiene kit

ref. F00935

~~€529,14~~

€465



Perio kit

ref. F00936

~~€670,97~~

€589



Implant Protect kit

ref. F02120

~~€735,54~~

€669



Endo Success Retreatment kit



ref. F00737

~~€702,36~~

€620



- Excellent view, full LED ring
- Avoid necrosis up to 120ml irrigation with peristaltic pump
- Improve durability of the handpiece, external irrigation
- Maximum hygiene, sterile single use or autoclavable lines
- Work don't play, use of big bottles (1l) possible
- Large choice of clinical application, 38 tips
- Tips II are compatible with IC2, Piezotome Cube, Piezotome2, Piezotome solo




piezotome cube

Dynamic power adaptation for
safe bone surgeries

~~€7.609,01~~
€5.650

Ref. F50100

Delivered with a multifunction footswitch, a Piezotome® CUBE LED handpiece, 1 handpiece holder, an Essential II kit (6 tips BS1S, BS4, SL1, SL2, SL3 et LC2), a Piezotome® dynamometric wrench, 1 irrigation pole, 2 autoclavable irrigation lines, 5 clips and 15 perforators.



Intraoral surgery ultrasonic control console.
Class IIa medical devices (GMED) - CE 0459
Manufacturer: SATELEC® a company of ACTEON® group - France

implant center 2

3 modes:

- Piezotome® (bone surgery).
- Newtron® (perio, endo,...).
- I-Surge™ (implant motor).

Ref. F27200 - Delivered with a LED motor, a LED motor cord, a progressive footswitch, 2 brackets, 2 handpiece holders, 2 sterilizable irrigation lines and 30 perforators.



~~€14.538,36~~
€9.599

The alliance of technologies
for powerful and secure
bone surgeries

Ultrasonic control console for conventional dentistry and intraoral dentistry and motor control console for dental implantology
Class IIa medical devices (GMED) - CE 0459 - Manufacturer: SATELEC® a company of ACTEON® group - France

Accessoires

Promotion valable only for ImplantCenter2

Newtron® LED handpiece



Ref. F12609
For Piezotome® 2 and ImplantCenter™ 2

Class IIa medical devices (GMED) - CE 0459
Manufacturer: SATELEC® a company of ACTEON® group - France

Newtron® LED cord



Ref. F12830
For Piezotome® 2 and ImplantCenter™ 2



~~€2.022,97~~
€1.489

Universal irrigation lines Sterile x10



~~€184,07~~
€145

as from 4 pieces

Ref. F59905 - Single-use with integrated cassette and flow control system.

Class Is medical device (GMED) - CE 0459
Manufacturer: SATELEC® a company of ACTEON® group - France

PIEZOTOME® LED handpiece + cord



Ref. F12815

Class IIa medical devices (GMED) - CE 0459
Manufacturer: SATELEC® a company of ACTEON® group - France

~~€1.635,16~~
€1.299



Bone Surgery™ II Kit*

Ref. F87509 - Tips BS1S II, BS2L II, BS2R II, BS4 II, BS5 II, BS6 II.

~~€1.432,35~~ **€1.090**



Sinus Lift II Kit*

Ref. F87519 - Tips SL1 II, SL2 II, SL3 II, SL4 II, SL5 II.

~~€1.077,92~~ **€990**



Universal II Kit*

Ref. F87539 - Tips BS6, LC1, LC2, CE1, CE3 et PZ1.

~~€1.393,47~~ **€1.090**



Essential II Kit*

Ref. F87528 - Tips BS1S II, BS4 II, SL1 II, SL2 II, SL3 II, LC2 II.

~~€1.432,35~~ **€1.090**



Intralift™ II Kit*

Ref. F87536 - Tips TKW1 II, TKW2 II, TKW3 II, TKW4 II, TKW5 II (Ø 3 mm).

~~€1.367,71~~ **€1.090**



Extraction™ II Kit*

Ref. F87546 - Tips LC1 II, LC2 II x2, LC2L II, LC2R II, Ninja® II.

~~€1.087,60~~ **€990**



Crest Splitting II Kit*

Ref. F87567 - Tips CS1 II, CS2 II, CS3 II, CS4 II, CS5 II, CS6 II.

~~€927,18~~ **€819**



Crown Extension II Kit*

Ref. F87554 - Tips BS6 II, CE1 II, CE2 II, CE3 II.

~~€737,70~~ **€619**

* Kits are delivered with an autoclavable metal tip holder and an autoclavable universal wrench.
Class IIa medical devices (GMED) - CE 0459 - Manufacturer: SATELEC® a company of ACTEON® group - France

* Kits are delivered with an autoclavable metal tip holder and an autoclavable universal wrench.

servotome

Delivered with a single footswitch,
a box with 10 electrodes,
a conductive bracelet.


Ref. F57210

~~€2.025,86~~
€1.699




Electrodes


Fulguration & coagulation





Coagulating incision



Excision



Incision



Patient's arm


	6 - 8	3 - 4
	6 - 8	3 - 4
	6 - 8	3 - 4
	5 - 7	1 - 3

SERVOTOME: Intraoral electrosurgery control console - Class IIb medical device (GMED) - CE 0459 - Manufacturer: SATELEC® a company of ACTEON® group - France
ELECTRODES: Class IIb medical device (GMED) - CE 0459 - Manufacturer: SATELEC® a company of ACTEON® group - France

air-n-go easy

All-in-one polisher for easy procedures

- AIR-N-GO® easy with KaVo* adapter Ref. F10100
- AIR-N-GO® easy with Sirona* adapter Ref. F10101
- AIR-N-GO® easy with W&H* adapter Ref. F10102
- AIR-N-GO® easy with Bien-Air* adapter Ref. F10103
- AIR-N-GO® easy with Midwest* adapter Ref. F10105



with a convertible AIR-N-GO® easy handpiece, a Supra nozzle, 1 CLASSIC powder bottle neutral (250 g), 1 AIR-N-GO® easy maintenance kit (replacement o’rings, syringe of silicone grease with brushes, cannula, cleaning drill).

*KaVo, Sirona, W&H, Bien-Air and Midwest are registered trademarks belonging respectively to the company of the same name.

air-n-go perio

air-n-go easy + air-n-go perio > ~~€2.081,16~~
€1.580



Kit PERIO easy
with a AIR-N-GO® easy nozzle and a bottle of
AIR-N-GO® PERIO glycine powder 160g.
Ref. F10120

Dental air-polishing handpiece - Class IIa medical device (GMED) - CE 0459.
Manufacturer: SATELEC® a company of ACTEON® group - France

AIR-N-GO® “CLASSIC” powder

4 bottles 250 g

- Peppermint Ref. F10250 ~~€105,48~~
- Lemon Ref. F10251 **€99.71**
- Raspberry Ref. F10252 each
- Cola Ref. F10253 **€74,90**
- Neutral Ref. F10254 as from
- Assortment Ref. F10255 4 pieces



AIR-N-GO® CLASSIC: Class IIa medical device - (GMED) - CE 0459 - Manufacturer: SATELEC® a company of ACTEON® group - France
AIR-N-GO® PERIO: Powder for subgingival treatment of periodontitis and peri-implantitis - Class IIa medical device - (GMED) - CE 0459
Manufacturer: SATELEC® a company of ACTEON® group - France

air-n-go perio



AIR-N-GO® “PERIO” powder
Ref. 10270
3 bottles 160 g

~~€130,39~~
€115

miniled standard

—€1.090,24—

€890

Ref. F02507

- **Fast**

Only 6 to 12 sec to polymerize
2mm composite.

- **Power**

1250mW/cm²
(Ø7,5mm light guide).

2000mW/cm²
(Ø5,5mm light guide - optional).

Highly Adaptable


The alliance of power & safety

- **Safety**

Base-charger includes a radiometer for light power
intensity check-up.

- **Optimal result**


3 modes to fit all clinical needs.



FAST MODE

Fast and easy


Full power during 10 sec
(bip at 5 sec) on regular
depth.



PULSE MODE

Ally rapidity & safety

10 flashes of 1 sec
(bip after 5 flashes)
with times off to
preserve dental pulp.



RAMPING MODE

No stress cracking

Progressive power
increase during 20 sec
(bip every 5 sec) for
a complete and safe
curing.



Simply Powerful

Outperforms in
power, efficiency and speed



~~€1.400,53~~

€1.090

Ref. F05214

miniled supercharged

• Visibility

Luminous screen with customizable modes and times, count curing time and load indicator.




• Power

2000mW/cm²
(Ø7,5mm light guide).

3000mW/cm²
(Ø5,5mm light guide - optional).

• Optimal result


3 modes to fit all clinical needs.



FAST MODE

The choice of rapidity


Full power during 3, 4, 5 or 10 sec, according to the composite thickness.



PULSE MODE

Limits temperature rise

5 or 10 flashes of 1 sec for a polymerization in security.



RAMPING MODE

Avoids composites retraction

Progressive power increase during 9 or 15 sec to ensure perfect sealing of restoration.

Dental curing light - Class I medical devices - CE
Manufacturer: SATELEC® a company of ACTEON® group - France

miniled ortho 2

~~€1.450,04~~

€1.090

Ref. F05220

• Customizable signal

An alert (beep, micro-flash, mixed or no alert) prompts the user to move to the next curing position every 4 sec.

• Power

3000mW/cm²
(Ø5,5mm light guide).

• Long lasting battery

MiniLED Ortho 2 can cure 15 patients at maximum power, no need to recharge it during the day.

• Ortho friendly

Customizable exposure times of 4, 8, 12 and 32 sec, specifically studied to move from one bracket to another or to cure a full arch.

• Efficacy & safety

Cure any type of material used orthodontics efficiently and quickly.

No noisy and unhygienic fan needed to cool the lamp.

Specialized Orthodontic Treatments

Cure in a simple gesture, in few seconds from
few brackets to a whole quadrant

Dental curing light - Class I medical devices - CE
Manufacturer: SATELEC® a company of ACTEON® group - France



Mouth mirrors

PURE REFLECT®

- Very bright image
- Resistant

For all purposes

RH

- A greyish reflection



TI

- A blueish reflection

FRONT SURFACE

Cross-sectional view






Image reflected directly on the surface of the glass, distortion-free.
NO DOUBLE IMAGE

BACK SURFACE

Cross-sectional view of an aluminium mirror



Reflecting surface under the glass. **DOUBLE IMAGE**



PURE REFLECT®

N°	Ref.	N.
4	MP3622PH	12
5	MP3624PH	12

~~€54,89~~
€39

Mirrors: Class I professional dental medical devices - CE
Manufacturer: PRODONT HOLLIGER ACTEON® Group - France

etik collagene

~~€90,45~~
€79

Hemostatic surgical dressing

Etik™ Collagene

25 unit-packaged sterile cylinders (in blister packs)
Ref. 281 970



Class III medical device (GMED) - CE 0459 - Manufacturer: PRODUITS DENTAIRE PIERRE ROLLAND - ACTEON® Group - France

softly™

Pre-measured capsules amalgamator

~~€496,85~~
€415

SOFTLY SAT 230V SCHUKO PLUG
Ref. 6.01.00112



Medical device classe I - CE - Manufacturer: SATELEC® Groupe ACTEON® (France)

Riskontrol® contributes to clinical protection



RISKONTROL® is also a GREEN solution

- Disposable tip
- Economic: simplify process to save time
- Recyclable plastic: environmentally friendly
- Cost-effective
- Flexible for a better use: respect of patient's anatomy
- Safety lock system: no risk of ejecting the tips




Sterilizable tips are not sufficient anymore

RISKONTROL® can be used for all procedures in any workflow




Easy to fit

1




It's installed.

2



It's fitted.

3



It's locked.

Discover our RISKONTROL® range

Adapt your tips to the taste
of your patients



by **art riskontrol** ~~€61,53~~
€32,50

- Ref. 201 700 - Mint
- Ref. 201 701 - Mandarin
- Ref. 201 702 - Liquorice
- Ref. 201 703 - Anise
- Ref. 201 704 - Black currant

Non-sterile disposable tips (x 250)

Protect entirely your air/water syringe
thanks to the plastic sheath



total protec
by **riskontrol**

~~€32,93~~
€29,50



- Ref. 201 750
- Ref. 201 751

Non-sterile disposable tips with plastic sheath (x 100)

**When you buy
10 boxes 1 adaptor FOC**

Add color
to your patients' care experience



~~€56,65~~
€32,50

by **classic riskontrol**

- Ref. 201 664 - White
- Ref. 201 665 - Blue
- Ref. 201 666 - Yellow
- Ref. 201 667 - Green
- Ref. 201 668 - Pink

Non-sterile disposable tips (x 250)

Choose a totally sterile
workflow



~~€38,13~~
€30,50

sterile
by **riskontrol**

Ref. 201 770


Sterile disposable tips (x 96)

Classic by RISKONTROL®, Art by RISKONTROL®, Total Protec by RISKONTROL®: Class I medical device. Sterile by RISKONTROL®: Class Is medical device. CE0459 (GMED). Read the instructions on the packaging carefully before use. All the compulsory information requested for a correct use of these devices is in the Product Feature Summary available on the laboratory's website. Manufacturer: PRODUITS DENTAIRE PIERRE ROLLAND - ACTEON® Group - FRANCE.

Use Expasyl™ Exact to improve patient comfort:

- 37 less pressure than traditional cords
- Preserved epithelial attachment
- Astringent effect - no bleeding
- Atraumatic sulcular opening
- No damage of the sulcus

(N) Force required for movement with retraction gels and compounds (in-vitro)



Adapted from American Association for Dental Research

Expasyl™ generates pressure that is 1.7 to 9.2 times greater than other retraction paste available on the market.

A TRIPLE ACTION



Kaolin clay



Aluminum chloride

1 ASTRINGENT EFFECT

helps to completely dry the sulcus with no blood or fluids.

2 OPEN THE SULCUS

effectively with no pain or damage to tissues thanks to the viscosity carefully calculated with our clinical expertise.

3 PRESERVE EPITHELIAL ATTACHMENT*

for a reduced risk of gingival recession and bone resorption.

*Canula must not be inserted into the sulcus

expasyl
exact



€144,33

€99

Paste for temporary sulcular opening
20 capsules (0,3g) Expasyl™ Exact - Ref. 261010



expasyl

€250,09

€139

Paste for temporary sulcular opening
20 capsules (1g) Expasyl™ - Ref. 261030



€111,70

€79

Straight - Ref. 261005
Bended - Ref. 261015

x-mind dc



The X-MIND® generators can accomodate any operatory configuration.

3 arms lengths are available : 41 cm (16.2''), 82.5 cm (32.5'') and 110 cm (43.5'') and can be:

- Top wall mounted
- Bottom wall mounted

Dental unit, ceiling and mobile version are also available.

With its adaptable wall plate, it is easy to replace an older generator with the X-MIND® whitout having to do additional drilling on the wall.

~~€4.611,95~~
€3.150

Available with Remote Wireless Control



No more wiring issues

Say goodbye to the complicated and expensive wiring work traditionnally needed to connect the remote button in your practice.

Wireless Button Kit X-MIND DC

Ref. 5.15.00135

All X-Mind DC devices delivered from April 2021 onwards

€249



Cost effective

Save on installation expenses with the easy-to-integrate wireless solution.



Enhanced convenience

Enjoy the freedom and flexibility of a wireless connection, making your practice more efficient and user-friendly.



Wall version	Top mounting	Top mounting square cone	Bottom mounting	Bottom mounting square cone
0,40 m extensor	Ref. 5.92.23137	Ref. 5.92.23167	Ref. 5.92.23130	Ref. 5.92.23160
0,80 m extensor	Ref. 5.92.23138	Ref. 5.92.23168	Ref. 5.92.23131	Ref. 5.92.23161
1,10 m extensor	Ref. 5.92.23139	Ref. 5.92.23169	Ref. 5.92.23132	Ref. 5.92.23162

All references for the X-MIND® DC

Choose from a variety of cones to better fit to your clinical needs



Long cone (31 cm - 12'')
Recommended for the parallel technique



Short cone (20 cm - 8'')
Recommended for bisector technique



Rectangular Long cone (31 cm - 12'')
Reduces the X-rays applied to the patient's skin by 50%.

Class IIb medical devices (IMQ) - CE0051 Manufacturer: de Götzen ACTEON® Group - Italy

x-mind unity

~~€5.767,06~~
€3.990

A sharp and contrasted image

- A contrasted and defined image
- Optimized protection for the patient
- The generator and sensor association

Fixing from above, arm 0.80 m, short cone.
Other configurations: see public price 2024.

Intraoral X-ray

- Very high quality
- Precise and flexible positioning
- Fast selection from traditional to digital mode
- Minimum exposure time: 0,02s
- Programmable timer
- Top or bottom mounting

Class IIb medical devices (IMQ) - CE0051 Manufacturer: de Götzen ACTEON® Group - Italy





sopix 2



A remarkable image quality, even for endodontics and the more complexes

SOPIX^{®2} provides a high resolution image by exploiting the latest technologies available on the market.



SOPIX^{®2} USB Sensor size 1
high definition with a sensor
holder, 10 sensor sheaths
Ref. S_802_0008

~~€6.920,49~~ **€3.990**

Also available in size 2
Ref. S_802_0009

~~€8.181,45~~ **€5.490**

Class IIa Medical devices (GMED) - CE0459
Manufacturer: Sopro ACTEON[®] Group- France

Free yourself from constraints
with the embedded system and its
detachable sensor part.



SOPIX^{®2} Plug-in kit S1

Ref. SK60-1000

~~€6.524,54~~ **€5.990**


SOPIX^{®2} Plug-in kit S2

Ref. SK60-1001

~~€8.580~~ **€7.490**

Size 2
31x42mm

Size 1
25x39mm



sopix 2 plug-in

The latest technology with no
compromise on image quality with
more safety, freedom and comfort.

Resolution: >18lp/mm

Class IIa Medical devices (GMED) - CE0459
Manufacturer: Sopro ACTEON[®] Group- France

x-mind optima 3D

The Smartest, The Smallest, The Ultimate.

In today's rapidly evolving field of dentistry, the integration of CBCT and Artificial Intelligence is transforming patient care, enabling dental professionals to make more accurate diagnoses and perform safer procedures.

Embracing this cutting-edge integration is not only elevating the standard of care but also optimizing efficiency and patient communication in modern dental practices.

X-MIND® optima 3D represents the perfect combination of advanced hardware technology, enclosed in a lightweight and compact design and next-generation software, RealGuide™ 5.4, bringing the benefits of Artificial Intelligence into every step of patient care, from a precise diagnosis to more effective treatments and patient communication.

Meet the X-MIND® optima 3D Family.

—€81.999—
€44.499

* Option:
The additional cost for
the 17x12 version is €3,000

X-MIND® optima 2D/3D wall 100-240V

2D/3D PAN device 5x5 to 12x10 wall installation
Ref. W1203005



—€93.499—
€54.999

* Option:
The additional cost for
the 17x12 version is €3,000

X-MIND® optima 2D/3D ceph 100-240V

2D/3D PAN CEPH device 5x5 to 12x10 floor installation
Ref. W1203010




—€83.499—
€45.999

* Option:
The additional cost for
the 17x12 version is €3,000

X-MIND® optima 2D/3D floor 100-240V

2D/3D PAN device 5x5 to 12x10 floor installation
Ref. W1203005.FL

New



Exceptional clarity.

At the heart of the X-MIND® optima 3D is **state-of-the-art detector technology**, which guarantees exceptional image clarity and outstanding anatomical coverage. Its **comprehensive range of Fields of View (FOV)**, from 5x5 to 17x12* cm (Ø x H), makes it suitable for a **wide range of dental procedures**, from general dentistry to orthodontics and maxillofacial surgery, while increasing tolerance for patient positioning and reducing the risk of marginal artifacts.



EXTRAORDINARY CLARITY

Get the most out of quality with the brand-new **ClearIQ** filter



OUTSTANDING VERSATILITY

From 5x5 to 17x12*, a tailored solution for every clinical need



TOP-LEVEL PERFORMANCE

Enhanced image depth with the new kV output range 60-90



MAXIMUM PATIENT COMFORT

Quick and easy acquisition workflow




AUTOMATIC POSITIONERS RECOGNITION

Maximum accuracy



ZERO FOOTPRINT


Maximum accessibility



* Option - The additional cost for the 17x12 version is €3,000

5x5 cm FOV

- Improves the **detection of root canals**, particularly in case of abnormal root canal shapes




8.5x5 cm FOV

- Useful for **implant surgery** and in **general dentistry** as a low dose exam as limited to a single jaw




8.5x9 cm FOV

- Useful for **complex implant surgery** as well as left or right TMJ diagnosis



12x10 cm FOV


- Suitable for most **oral surgery cases** including multiple implant surgery



17x12* cm FOV

New

- Represents the perfect solution for **complex orthognathic cases**, **zygomatic implant planning** and **orthodontic surgery**



* Option - The additional cost for the 17x12 version is €3,000



AI
for life



ACTEON CLUB
Welcome leaflet + Webkey



Workstation and monitor



acteon
imaging suite
v 5.4



x-mind
optima 3D



RealGUIDE™
UNIVERSAL OPEN SYSTEM
v 5.4 PRO

AIS 5.4

RealGuide™ 5.4 PRO

X-MIND® optima 3D options offer

Option	Ref.
RealGuide™ 5.4 DESIGN (1 year license)	W1000026
17x12* FOV	W1203107

* Option - The additional cost for the 17x12 version is €3,000

x-mind
prime 2


~~€18.395,52~~
€15.999

5
YEAR
WARRANTY**

X-MIND® Prime2 2D wall mounted

The ACTEON®'s smallest extraoral 2D imaging solution ever for the **maximum accessibility in any context.**

Ref. W1200001 - W1200001.PB*



PANORAMICA STANDARD HD - 80 kV 9.0 mA 14.0 s - 157.2 mGy cm2



~~€54.195,29~~
€29.999

5
YEAR
WARRANTY**

X-MIND® Prime2 2D ceph

Cephalometric add-on for a full 2D system still encompassing the same philosophy of prime, **ideal for pediatric applications.** Included 2 sensors.

Ref. W1200009.DS - W1200009.BP*



LATERO-LATERALE 30 X 24 CM HD - 76 kV 8.0 mA 15.1 s - 34.8 mGy cm2

Class IIb medical devices (IMQ) - CE0051 Manufacturer: de Götzen ACTEON® Group - Italy
* ONLY for Belgium with lead
**2 year full factory warranty + 3 year for Tube and Generator only!

PSPiX²

is AUTONOMOUS


The fully automated process free yourself from constraints

Just drop the plate inside the slit, the PSPiX² will automatically detect the size, scan the plate, optimize the image, and send it instantly to your imaging software.

In the meantime, it will automatically erase and eject the plate for an immediate reuse or switch to standby mode.

Scan
simply &
quickly

Fully automated process



—€8.358,04—

€ 6.490

Ref. S_702_0001

1

Auto-access:

The door only opens when the imaging plate is inserted.

2

Auto-detect:

PSPiX² automatically detects the size of the inserted imaging plate.

3

Auto-scan:

PSPiX² scans and optimizes the image before displaying it into the right workstation.

4

Auto-eject:

The imaging plate is automatically erased and ejected for an immediate reuse.

5

Eco-mode:





PSPiX² automatically switches to economic mode when not in use for an extended period of time.

Protective bag for ERLM*

For **PSPIX®2** and **PSPIX®**
(previous generation)

~~€84,14~~
€57,40



	PSPIX® (previous generation)	
	Box of 250	
Size 0		700338L
Size 1		700339L
Size 2		700340L
Size 3		700341L





Protections for ERLM*

Protective cover for
PSPIX® and **PSPIX®2**
(previous generation)

~~€47,12~~
€26,25



Compatible with	PSPIX²	PSPIX
	✓	✓
PSPIX®	✓	✓
SOREDEX® Digora	✓	✓
INSTRUMENTARIUM® Express	✓	✓
DÜRR DENTAL® Vistascan®	✓	✓
KAVO®	✓	✓





	PSPIX® (previous generation)	
	Box of 250	
Size 0		700342L
Size 1		700343L
Size 2		700344L
Size 3		700345L

Capteurs ERLM*

For **PSPIX®2**

~~€131,80~~
€95







	PSPIX®2	
	Box of 2	
Size 0 (22 x 35mm)		990215
Size 1 (24 x 40mm)		990216
Size 2 (31 x 41mm)		990217
Size 3 (27 x 54mm)		990218

Protective bag for ERLM*

For **PSPIX®2**

~~€98,88~~
€70



	PSPIX®2	
	Box of 300	
Size 0		700546
Size 1		700547
Size 2		700548
Size 3		700549

- Protection against scratches, light and cross contamination
- Ergonomic and easy to use
- Economic solution

* Luminescent Radio Memory Display for digital phosphor plate systems

PSPIX® imaging plate - Class IIa Medical devices (GMED) - CE0459
Manufacturer: Sopro ACTEON® Group - France
Bag & cover: Class I medical devices - CE
Manufacturer: Sopro ACTEON® Group - France



PSPiX²

is ADAPTABLE


It fit into any practice and designed for everyone

Compact Size

Created with cutting edge technology and design, the PSPiX² is the smallest and most compact PSP scanner on the market.

Optimal Protection

The autoclavable option allow removable part to comply with thermal washer-disinfectors, the latest disinfection standards.



So Simple

The large colour touchscreen provides intuitive indications for quick and easy use.



AVAILABLE



SCANNING



BOOKED

Use it as your personal assistant...

ACTEON[®] has created the first personal imaging plate scanner on the market:

Get your images immediately to experience high productivity.

With PSPiX² chairside, you significantly save time:


- You are the single-user of the scanner
- You do not need to move to another room to scan your images
- You get your images in a few seconds
- You are therefore more efficient and more dedicated to your patients' care

... or in collaborative mode

PSPiX² has also been designed to be shared by up to 10 workstations

Save time with the new “Click and Scan” concept of PSPiX², the simplest way to gain efficiency.


Select your workstation on the touchscreen, insert your imaging plate and let PSPiX² do the rest...



The thin and flexible imaging plates, available in 4 sizes, will allow you to easily scan all morphologies...




... and all clinical cases




C50

Daylight



Transcript of the **natural colors** for a better interpretation of the dental tint and helps to identify pathologies and with post-treatment control.

Daylight*



High contrasts emphasizing anatomic fissures and details of the dentine structure / depth.

- Full HD: 1920x1080
- Ergonomic design with small tip (easy to use with small mouth/children)
- Autofocus
- Diagnostic feature: dental plaque and caries detection (C50 Fluo)

€3.177,85

€2.699

(daylight)

Ref. S_C50_0002 DAY LIGHT

€4.574,56

€3.499


(fluorescence technique)

Ref. S_C50_0001 FLUO




4 Modes

Daylight




Transcript of the **natural colors** for a better interpretation of the dental tint and helps to identify pathologies and with post-treatment control.

Daylight*




High contrasts emphasizing anatomic fissures and details of the dentine structure / depth.

Perio



Enhancement of **dental plaque and gingival inflammation** thanks to chromatic amplification.

Cario



Caries **detection by fluorescence technique** by LEDs.

5 Visions


Macro



Intraoral



Portrait



Tooth



Smile



x-mind optima3D

EASY-JOIN QR Code

SCAN ME!



 acteon club